

Integrated, closed-loop
sampling with aScope™

Ambu® aScope BronchoSampler™

Streamline procedures and
protect samples with aScope
BronchoSampler, a fully-
integrated, closed-loop solution
only for aScope™ 4 Broncho.

A simplified sampling process

A one-of-a-kind, purpose-built system, the aScope BronchoSampler improves workflow efficiency by allowing a single clinician to easily perform bronchoscopic sampling as part of BAL or BW procedures – from system assembly to sample procurement.

aScope BronchoSampler also features a vacuum bypass that further simplifies procedures by enabling you to change between suction and sampling without a suction tube switch.

Fewer lost and contaminated samples

As a single-use, sterile, closed-loop system aScope BronchoSampler reduces the risk of sample loss and guarantees sample quality from start to finish. By minimising the risk of sample contamination or exposure, aScope BronchoSampler also supports accurate and timely diagnostics as well as the safety of both patients and clinicians.

Innovative, integrated system

aScope BronchoSampler is designed to integrate seamlessly with your aScope 4 Broncho. It comes with all the critical components that you need, allowing you to effectively and efficiently perform BAL and BW procedures. And because it's always ready when you are, aScope BronchoSampler minimises the significant risks to patients that can result from equipment and workflow delays.



aScope BronchoSampler™


– everything you need in a sterile, portable and pre-packed tray

Lid closes to prevent sample contamination

Flip top to easily disengage sealing connection and remove sample container

Flow switch for easily changing between suction and sampling

Interface for sealing connection



Transparent container allows for monitoring of sample volume and texture

Attachment lock to aScope 4 Broncho

Suction connection

Ordering specifications

Item no	Quantity
480200000 aScope bronchoSampler™, specimen sampling system	5
480300000 aScope bronchoSampler™ SC, sample container	10

Ambu Ltd.
8 Burrell Road
St Ives
Cambridgeshire
PE27 3LE
Tel 01480 498403
Fax 01480 498405
uksales@ambu.com

For more information
Visit ambu.co.uk

4496480100 - V01 - 2018/09 -
Ambu A/S. Technical data may be
modified without further notice.